Vacotron 460 / 560

Bedieningshandleiding

Operating Instructions

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Mode d'emploi

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Istruzioni d'uso

Copyright:

**Enraf
Nonius**

Enraf-Nonius B.V.

P.O. Box 12080

3004 GB Rotterdam

The Netherlands

Tel.: +31(0)10 - 20 30 600




Fax: +31(0)10 - 20 30 699

Order number 1490.750-45

June 1999

Vacotron 460 / 560

**Bedieningshandleiding
Operating Instructions
Gebrauchsanweisung
Mode d'emploi
Instrucciones de manejo**

Vacotron 460

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TABLE OF CONTENTS

Foreword	13
General	13
Constant or pulsed vacuum	13
Electrodes	13
Smart connections for simple operation	13
Finally	13
Introduction	14
General	14
Installation	14
Connection	14
Electrolytic effects	14
Current density	14
Electromagnetic interference	14
Technical maintenance	14
Product liability	15
International safety standard IEC 601-1	15
Installation	15
Attaching the water reservoir (rear side)	15
Connecting the electrotherapy unit (rear side)	15
Connecting the vacuum hoses (front side)	15
Patient cables (front side)	15
Mains cable	15
Controls	16
Operating sequence	17
VACUUM ELEKTRODES	17
Switching on	17
Positioning the electrodes	17
Correcting the vacuum	17
Selecting the current	17
Ending treatment	17
Flexible rubber electrodes	17
Contra-Indications	17
Maintenance by the User	18
Treatment unit	18
Electrodes and sponges	18
Electrode cables	18
Water reservoir	18
Fault conditions	18
Lamps fail to light	18
Changing a fuse	18
Water in the pump	18
No patient current	18
Patient current too low	18
Specifications and Technical data	19
Ordering data	19
Standard accessories	19
Additional accessories	19

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General

The Vacotron 460 / 560 pulsed vacuum unit provides fast and simple electrode positioning using vacuum electrodes. Furthermore, it is extremely simple to switch between treatment with vacuum electrodes and treatment using flexible rubber electrodes. As a result, the Vacotron 460 / 560 is an ideal aid to electrotherapy.

When vacuum electrodes are applied, effective use is made of the whole electrode area, as the electrodes make good contact with the skin. The massage effect resulting from the pulsed vacuum ensures a good flow of blood through the skin under the electrode. This reduces the resistance of the skin, increases the effectiveness of the stimulating current, and reduces the treatment time.

Constant or pulsed vacuum

When a constant vacuum strength is used, the vacuum pulse switches off as soon as the selected vacuum strength is reached, making the Vacotron 460 / 560 virtually silent in operation. A control system ensures that the vacuum remains constant. The pump switches on for short intervals during the treatment to maintain the vacuum.

When pulsed vacuum is used, the pump is switched on and off intermittently. The pulsating rhythm of the vacuum massage is continuously or pulsating (two possibilities) adjustable. Pulsating vacuum is very useful for patients who are sensitive to the electric current, as the sensation of the current tends to be suppressed.

Electrodes

There is a choice of large and small electrodes. The areas of the electrodes correspond to those of the 4 x 6 cm and 6 x 8 cm flexible rubber electrodes. The vacuum electrodes are sufficiently flexible to ensure optimum contact with the skin, but rigid enough to prevent any changes in the contour of the part being treated, allowing full advantage to be taken of the massage effect of the pulsed vacuum.

Smart connections for simple operation

The Vacotron 460 / 560 has been designed to ensure easy operation. This is reflected in the design of the electrode connections. The Vacotron 460 / 560 has connections for both vacuum electrodes and flexible rubber electrodes, providing a permanently available choice between the two types without the need to change the cables. The vacuum electrodes are automatically selected as soon as the vacuum is set.

Finally

The Vacotron 460 / 560 has been designed and manufactured to ensure the highest level of safety, and the unit fully complies with the stringent IEC 601-1 standards. We are confident that your unit will continue to give satisfaction over many years of use. Nevertheless, if you have any queries or suggestions, please contact your authorized distributor.

General

It is important that you read these operating instructions carefully before using the Vacotron 460 / 560. Please make sure that these instructions are available to all personnel who operate the equipment.

Pay attention to the following before using the Vacotron 460 / 560:

1. Keep yourself informed of the contra-indications (see page 17).
2. The apparatus may not be used in close proximity (i.e. less than 2 metres) to shortwave or microwave equipment.

The manufacturer cannot be held responsible for the results of using this apparatus for any purposes other than described in these operating instructions.

Installation

- Do not install the unit in a location near a heat source such as a radiator.
- Avoid exposure to direct sunlight, rain, excessive dust, dampness, mechanical vibrations and shocks.
- This unit may not be used in so-called "wet rooms" (hydrotherapy rooms).

Should any liquid or foreign object get into the housing of this equipment, turn off the power switch immediately, unplug the unit from the wall socket and have it checked by an authorized expert.

Connection

- Connection to the electrical mains supply must comply with the national requirements regarding medical rooms.
- Before connection of this apparatus to the supply, check that the mains voltage and frequency stated on the type plate of this apparatus correspond to that of the mains supply.
220-240 V model for countries with mains voltage of 220, 230 or 240 volt.
115 V model for countries with mains voltage of 110 or 120 volt.
- This apparatus may only be connected to an earthed wall socket.
- The unit should only be used with the mains cable supplied with this equipment.
- A potential equalisation cable may be connected to this unit, if prescribed in the national regulations.

Technical maintenance

We recommend to have the unit checked annually. This may be done by your supplier, or by an other agency, authorized by the manufacturer. It is also recommended that a record of the service history is kept for all activities relating to service and maintenance. In some countries this is even obligatory.

No attempt should be made to open the unit. Maintenance and all repairs should be carried out by an authorized agency.

The manufacturer will not be held responsible for the results of maintenance or repairs by unauthorized persons.

Electrolytic effects

Electrolysis occurs under the electrodes when Galvanic or currents with a DC component are applied. Because the greatest concentration of electrolytic by-products caused by ion migration occur under the electrodes, we recommend the use of the supplied sponges to keep the effects to a minimum. Make sure that the sponges are kept well moistened.

Current density

The particular standard for Therapeutic Electrical Stimulation Equipment is IEC 601-2-10 (BS5724 part 2.10). The maximum permissible Current Density is 2 mA R.M.S. per cm².

To verify that this value is not exceeded in electrotherapy, divide the current output in mA by the effective electrode area in cm²; a value of less than two should result.

Note that placing the electrodes should be done carefully and solidly, and that the use of small electrodes in combination with high intensities may cause skin irritations or burns.

Electromagnetic interference

Simultaneous connection of a patient to High Frequency surgical equipment and an electrical stimulator must not be done as this may result in burns at the site of the stimulator electrodes.

Operation in close proximity (i.e. less than 2 metres) to working shortwave or microwave therapy unit may produce instability in the Vacotron 460 / 560 output.

To prevent Electromagnetic Interference, we strongly advise using separate mains groups (phases) for the Vacotron 460 / 560 and the short-wave or microwave equipment, and a distance of at least 2 metres is kept between the Vacotron 460 / 560 and the shortwave or micro-wave equipment. Ensure that the mains cable of the shortwave/microwave unit does not come near the Vacotron 460 / 560 or the patient.

If problems with Electromagnetic Interference persist, please contact your authorized distributor.

Product liability

A law on Product Liability has become effective in many countries. This Product Liability law implies, amongst other things, that once a period of 10 years after a product has been brought into circulation has elapsed, the manufacturer can no longer be held responsible for possible shortages of the product.

International safety standard IEC 601-1

This equipment meets the requirements of the International Standard for the safety of medical electrical equipment IEC 601-1. A list of countries along with their own national codes that have incorporated IEC 601-1 into their national standards:

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(see fold-out sheet at front)

Attaching the water reservoir (rear side)

Attach the water reservoir. This can be made easier if a thin layer of Vaseline is applied to the connection 4.

The reservoir also serves as a compressor chamber, and must be free from air leaks.

Connecting the electrotherapy unit (rear side)

The Vacotron 460 / 560 is suitable for use with single-channel or two-channel electrotherapy units. The unit may be of any type or make, but must meet the IEC 601-1 standard for type BF equipment. All ENRAF-NONIUS electrotherapy units comply with this standard.

Connect the electrotherapy unit to the connection sockets 2 and 3 at the rear of the Vacotron 460 / 560, using the patient cables supplied with the electrotherapy unit. The connections are placed next to each other, and are identified as channel 1 and channel 2.

For currents with a DC component (e.g. diadynamic currents) connect the positive pole (red) to the left hand socket and the negative pole (black) to the right hand socket. For alternating currents, connections are at choice. Consult the Operating Instructions of the electrotherapy unit for this information.

Connecting the vacuum hoses (front side) channel 1

- Connect the red hose to the left hand socket 8 (+).
- Connect the black hose to the right hand socket 8 (-).

channel 2

- Connect the red hose to the left hand socket 10 (+).
- Connect the black hose to the right hand socket 10 (-).

Patient cables (front side)

Connect the patient cables supplied to connection 15 for channel 1 and connection 17 for channel 2.

Mains cable (left hand side)

Connect the mains cable to connection 7. The Vacotron 460 / 560 is now ready for use.

1 Indication plate

Provides information on the unit, such as type and serial number, as well as connection data such as mains voltage and current consumption.

2 Connection 1 for electrotherapy unit

This connection is used for connecting one channel of the electrotherapy unit to the Vacotron 460 / 560. Use the patient cable supplied with the electrotherapy unit (see section on Installation).

3 Connection 2 for electrotherapy unit

This connection is used for connecting the second channel of the electrotherapy unit to the Vacotron 460 / 560. Use the patient cable supplied with the electrotherapy unit (see section on Installation).

Only type BF equipment may be connected to the Vacotron 460 / 560. The very small leakage current of this type of equipment ensures absolute safety.

All ENRAF-NONIUS electrotherapy units are of type BF.

4 Connection for water reservoir

This connection is used for connecting the water reservoir. No vacuum can be created unless the water reservoir is connected. The water reservoir ensures that any excess water from the sponges is collected, and cannot be blown out through the pump. The reservoir should be emptied daily. Please refer to the section on Maintenance.

5 Air outlet

The air outlet is provided with a sound-absorbent filter. Replace the filter once a year. Please refer to the section on Maintenance.

6 Potential equalization point

For use in areas where local regulations require connection of a potential equalization cable. An earth (ground) cable and plug can be supplied.

7 Mains switch, mains input and fuses

The mains switch is used to switch the Vacotron 460 / 560 on and off. The mains cable is connected to the mains input, which is provided with two fuses. The Vacotron 460 / 560 need only be switched on if vacuum electrodes are used. If the unit is switched off, the flexible rubber electrodes are selected automatically.

8 Vacuum channel 1

Left hand connection: the red vacuum hose (+).
Right hand connection: the black vacuum hose (-).
Vacuum channel 1 is connected to one channel of the electrotherapy unit via connection 2.

9 Vacuum electrodes indicator lamp

Lamp is lit when the vacuum electrodes are operative.

10 Vacuum channel 2

Left hand connection: the red vacuum hose (+).
Right hand connection: the black vacuum hose (-).

Vacuum channel 2 is connected to the second channel of the electrotherapy unit via connection 3.

Connection of any cables other than the vacuum cables can adversely affect the safety of the patient and the functioning of the unit, and is therefore not permitted.

11 Vacuum regulator

This regulator sets the vacuum level, and is also used to switch the vacuum system on and off.

No vacuum:

Turn the regulator fully to the left.

The vacuum electrodes are inoperative. Indicator lamp 9 is out. The patient cables connected to connections 15 and 17 are operative.

Vacuum:

Turn the regulator to the right. The indicator lamp 9 by the vacuum electrodes is on, indicating that the vacuum electrodes are now the operative electrodes. Further rotation to the right increases the vacuum.

12 Key for constant vacuum

After pressing this key the Vacotron will supply constant vacuum (an indicator lamp lights up).

13 Key for pulsed vacuum 0,5:0,5

After pressing this key the Vacotron will supply pulsed vacuum with a fixed pulse:pause ratio of 0,5:0,5 seconds.

14 Key for pulsed vacuum 1:1

After pressing this key the Vacotron will supply pulsed vacuum with a fixed pulse:pause ratio of 1:1 second.

15 Patient cable connection channel 1

Connect one of the patient cables supplied for use with the flexible rubber electrodes to this connection. This connection is selected when no vacuum has been set. Channel 1 is connected to one channel of the electro-therapy unit via connection 2.

16 Flexible electrodes indicator lamp

Lamp is lit when the flexible rubber electrodes are operative.

17 Patient cable connection channel 2

Connect the second patient cable supplied for use with the flexible rubber electrodes to this connection. This connection is selected when no vacuum has been set. Channel 2 is connected to the second channel of the electrotherapy unit via connection 3.

Connection of any cables other than the cables for the flexible rubber electrodes can adversely affect the safety of the patient and the functioning of the unit, and is therefore not permitted.

OPERATING SEQUENCE

VACUUM ELECTRODES

Switching on

- Switch on the Vacotron 460 / 560 and the electrotherapy unit.

Positioning the electrodes

- Choose the appropriate electrode size for the body part to be treated.
- Dampen the sponges. Squeeze them slightly to remove excess water, and place them in the vacuum electrodes.
- Select constant vacuum, this will ease the positioning of the electrodes.
- Place the vacuum electrodes on a smooth surface (e.g. on a table).
- Turn the vacuum intensity up a little, so as to create a vacuum in the electrodes.
- Now move the electrodes one by one from the smooth surface to the body of the patient.
- Set the vacuum intensity to the desired value.
- Set the treatment time.

The treatment starts.

Correcting the vacuum

- Correct the vacuum if necessary with regulator 11. A low vacuum is sufficient.
- If required, select the pulsed vacuum with the keys 13 and 14.

Selecting the current

- Select the required current type, and gradually increase the current. Please refer to the Operating Instructions supplied with the electrotherapy unit.

Ending treatment

- Reduce the current of the electrotherapy unit to zero.
- Turn the vacuum regulator back to position 2, and remove the vacuum electrodes from the skin.
- Turn the vacuum regulator back to zero.

FLEXIBLE RUBBER ELECTRODES

- When using the flexible rubber electrodes, the Vacotron 460 / 560 need not be switched on.
- If the Vacotron 460 / 560 is switched on, the vacuum regulator should be in the zero position. The indicator lamp 16 near to the patient cables connections will light up.
- Treatment continues as described in the Operating Instructions of the Electrotherapy unit.

N.B.

During treatment with the flexible rubber electrodes, simultaneous use of the vacuum electrodes is not possible. The vacuum regulator should remain in the zero position during treatment with the flexible rubber electrodes.

CONTRA-INDICATIONS

Supplementary pulsed vacuum therapy can be applied in almost every situation where low-frequency electrotherapy is indicated.

The use of the vacuum electrodes is only contra-indicated in the case of inflammatory processes, or where the risk of bleeding at the electrode sites cannot be excluded.

Treatment unit

First switch off the unit and disconnect the mains cable. The unit can be cleaned with a damp cloth. Use lukewarm water and, if required, a liquid household detergent. Do not use abrasive cleaners.

Electrodes and sponges

The vacuum electrodes and sponges should be cleaned with lukewarm water. In the case of persistent dirt, and for disinfection, a 70% alcohol solution may be used.

Sponges should be replaced regularly. It is recommended to keep sponges and a spare electrode in stock.

Calcium scale can be deposited on the metal surfaces of the electrodes. This has an insulating effect. In order to maintain optimum conductivity, these surfaces should be regularly cleaned and polished.

Electrode cables

Check the cables regularly for damage and/or poor contacts. It is recommended to keep a spare electrode cable in stock.

Water reservoir

The Vacotron 460 / 560 is provided with a water reservoir, designed to prevent water from dripping into the pump. The reservoir also serves to collect excess water from the sponges. The water reservoir should be checked after each treatment, and emptied and cleaned if necessary. The water reservoir should be emptied at least daily. If calcium scale appears, it can be removed with a domestic descaling agent.

The following registered products may be used for disinfecting the water reservoir:

BAKTOLAN to 5%
CHINOSOL to 1%
CHLORAMIN solution
ELMOCID Gamma to 2%
MEFAROL to 1%
MERCKOJOD to 1%
MERFEN
PERHYDROL
PERODIN
SAGROTAN to 2%
ZEPHIROL to 5%.

Lamps fail to light

Check the mains cable and the mains fuses.

Changing a fuse

There are two fuses in the mains input. Remove the mains cable and open the fuse holder. Pull out the fuse holder and replace the defective fuse with one of the same type (T1.6H250V). Push the fuse holder back in and reconnect the mains cable.

Water in the pump

Water in the pump can be recognized by a 'spluttering' sound. If water gets into the pump (e.g. because the reservoir has not been emptied on time) proceed as follows:

- Remove the reservoir and empty it;
- Remove the air outlet filter 6;
- Run the Vacotron 460 / 560 for about five minutes without the reservoir and filter (constant, maximum vacuum);
- The pump will then blow the water drops out of the system;
- Replace filter and empty reservoir.

No patient current

Check that the correct electrodes have been selected. If the vacuum regulator is in the zero position, the flexible rubber electrodes are selected automatically. As soon as a vacuum is set, the vacuum electrodes are selected.

Check operation of the electrotherapy unit using the operating instructions provided (if necessary, with the Vacotron 460 / 560 disconnected).

Check the cables from the electrotherapy unit to the Vacotron 460 / 560. The connections for channels 1 and 2 may have been inadvertently exchanged.

Patient current too low

The electrical resistance of the electrodes is too high. This can be due to the sponges not being damp enough, or to a layer of calcium scale or corrosion on the electrodes.

Vacuum	: 0-600 mbar, continuously adjustable
Pulsed vacuum:pulse	: pause = 0,5:0,5 seconds, 1:1 second
No. of vacuum electrodes	: max. 4
Mains supply	: 220-240 V~ ± 10%, 50/60 Hz 115 V~ ± 15%, 50/60 Hz
Current consumption	: max. 0.13 A (at 230V) max. 0.26 A (at 115V),
Mains fuses	: 2 x T1.6H250V
Patient leakage current	: < 10 µA (IEC requirement ≤ 100 µA)
Ditto, first fault condition	: < 10 µA (IEC requirement ≤ 500 µA)
Earth (ground) leakage current	: typically 30 µA (IEC requirement ≤ 500 µA)
Ditto, first fault condition	: typically 50 µA (IEC requirement ≤ 1000 µA)
Safety class	: I* type BF**, according to IEC 601-1
Dimensions	: 35 x 36 x 11 cm
Weight	: 5.7 kg
Safety	: TÜV Rheinland This device is in conformity with the European Directive (93/42/EEC).

*I The equipment has a safety earth (ground) connection, and must be connected to an earthed (grounded) wall socket.

**BF: The equipment has a floating patient circuit. The Vacotron 460 / 560 meets the requirements of IEC 601-1 and IEC 601-2-10.

Copies of the test reports are available on request.
Technical modifications reserved.

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